95060001 CA RTM 2018-11 18-63

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Instructions for Use Revision and Tumour System

The present Instructions for Use apply to the REVISIO® RTM revision and tumour system.

Appendix 1 contains an implant overview including all relevant item numbers“.

Improper implantation of the prosthesis can lead to failure, implant failure, patient injury, bone loss, or unsatisfactory aesthetic results. The products may only be used by trained physicians and surgeons who are familiar with the implant system. Use of these products requires specialised knowledge and skills in endoprosthesis/surgical implantology. AQ Implants GmbH offers specific training and technical consultation on the use of its products and provides detailed information for appropriate implant selection as well as preoperative planning and surgical technique. Our products are sold only to physicians and healthcare professionals or on their behalf. Technical consultation on our products will be performed verbally or in writing, by electronic media, and/or by demonstration. Consultation reflects the latest state of science and technology at the time of the product launch. It does not relieve the user from the obligation to personally verify suitability of the product in terms of intended purpose, indications, and procedures. Handling and implantation of the product are beyond our control and are the responsibility of the user. Any liability for damages shall be excluded.

1. Product description and intended use

The REVISIO® RTM endoprosthesis system was developed for complex applications in tumour surgery. Due to the high modularity of the REVISIO® RTM system, the individual components of the proximal and distal femur as well as the diaphysis up to the complete femur can be replaced. The REVISIO® RTM implant system also features an arthrodesis spacer for joint arthrodesis of the knee joint as well as the REVISIO® RTM condylar replacement to maintain the mobility of the knee joint. The REVISIO® RTM endofemoral module which is introduced through the femur connecting hip and knee implant represents an expansion of the REVISIO® RTM system. It was developed to minimise bone loss in the femur by preserving the diaphyseal femur. The modular implant components are available in various diameters and lengths ensuring an adequate reconstruction of the individual defect. The individual modules are connected by means of tapered clamping and additional fixation via a specific screw system.

The combination of implant components depends on the indication and is the responsibility of the surgeon.

2. Handling

Keep the implant in the unopened original packaging until it is used. When unpacking the implant, please ensure that the protective packaging, i.e. cardboard and peel film, is unopened and intact. Any damage may impair the sterility of the implant. In the event of a broken seal or damaged cardboard and peel film, it must be assumed that the contents are no longer sterile.

Caution:

If the seal is broken, the product cannot be returned.

When unpacking the implants, please verify the details on the packaging against the package contents (item number and size) as well as expiry date. Do not implant components after the expiry date stated on the package.

Caution:

Make sure that all aseptic standards are observed during removal of sterile implant components.

Take appropriate measures to avoid that the implant comes into contact with objects or substances that might damage its surface. This applies particularly for coated implants (calcium phosphate [CaP], Titanium-Niobium [TiNb] etc.) e.g. for the contact with body fluids (sweat, blood, synovia, etc.). For ceramic implants this also applies for the direct contact with skin. Please note the enclosed user information and surgical techniques. Visually inspect implant and coating (where applicable) for damage prior to implantation. Damaged implants or implant components must not be implanted. Protective caps or similar devices for the protection of sensitive implant areas are removed immediately before use. Implants must not be processed mechanically, altered or modified, unless this is explicitly required by the surgical technique or the manufacturer has given his written consent.

Caution:

Products that have been implanted must not be reused. The same applies for implants that came into contact with foreign body fluids.

Removal, handling and return of revised implants for scientific evaluation must be performed in accordance with ISO 12891-1. Please observe the general information of ISO 8828 on the care and handling of orthopaedic implants! For specific items, please observe the manufacturer's instructions (handling, sterilisation, compatibility, etc.).

3. Packaging and sterility

All implants are packed individually in sterile protective packages. Sterile implants have been sterilised with a minimum gamma radiation dose of 25 kGy.

Caution:

Any implants labelled „sterile“ are only sterile, if the package has not been opened.

Caution:

Implants must not be resterilised. Any reprocessing of non-implanted components with opened packages or expired sterile date must only be performed by the manufacturer, as several validated processes have to be performed again.

These instruments can be autoclaved using standard procedures (autoclave according to DIN EN 285, procedure according to DIN EN ISO 17665-1 and DIN EN ISO 14937), provided the required sterilisation temperatures and cycle times are observed. Prior to resterilisation, all instrument components should be disassembled as far as possible. Please refer to the „Reprocessing guidelines for the reprocessing of reusable surgical instruments“ of AQ Implants GmbH for detailed information.

Note:

Implants for use with or without bone cement are marked with the following abbreviations:

- Implants intended for use with bone cement: C (cemented)
- Implants intended for use without bone cement: CL (cementless)

4. Surgical planning

Preoperative planning is based on the radiographic evaluation. Appropriate X-ray templates are available. Moreover, all products of AQ Implants GmbH are available in the digital planning system mediCAD®.

For the preparation of the host bone and adjustment, repositioning as well as implantation of the implants, specific instruments and trial implants corresponding to the dimensions of the implants are available. Please don't use other trial instruments or implants.

The appropriate surgical technique is selected in accordance with the latest state of science and technology. In addition, AQ Implants GmbH provides surgical techniques upon request which are constantly revised in view of new findings and understanding.

In order to achieve optimal primary fixation as a prerequisite for long-term stability of the prosthesis, the following points should be observed:

- Preservation of adequate bone structure / wall thickness during the preparation of the host bone
- Selection of the appropriate implant size
- Adequate cleaning of the host bone prior to implantation
- Application of moderate force during implantation to avoid bone fissures or implant damage.
- Compatibility of the individual components
- Observance of surgical technique, user notes, and warnings
- Clean and dry surfaces of tapered connections
- Firm coupling of tapers achieved by using insertion and impaction instruments
- Avoid damage to the components
- Avoid foreign particles, such as bone or tissue debris or cement particles between functional surfaces (e.g. tapers, holes, articular surfaces)
- Use of locking screws in screw-retained components (e.g. fixation of the iliac peg)

In addition, the following points should be observed for cemented All-Poly Cups:

- Observance of the processing specifications for bone cement (PMMA), mixing ratio as well as polymerisation temperature and time in accordance with the manufacturer's instructions
- Appropriate mixing of polymer and monomer
- Application of bone cement of appropriate viscosity without blood and air inclusion
- Uniform cement mantle around the implant
- Appropriate press-fit and resting of the implant during polymerisation
- Thorough removal of excessive bone cement from bony edges
- Thorough rinsing of the surgical site prior to reposition and wound closure to avoid cement particles and bone fragments from getting trapped between joint surfaces or functional and clamping surfaces, holes, etc.

Caution:

The cement mantle should fully enclose the implant in the fixation area. An incomplete cement mantle can lead to early loosening of the implant.

Note:

Bearing surfaces and tapers of the prostheses can be damaged by surgical instruments or bone cement particles. In order to guarantee the safe and reliable function, they have to remain intact. Therefore, these areas should be covered with surgical drapes during implantation.

Note:

The tapers of the prostheses can be damaged by surgical instruments or bone cement particles. In order to guarantee the safe and reliable function, they have to remain intact. Therefore, these areas should be covered with the supplied taper covers and surgical drapes during implantation.

Bearing surfaces of the prostheses can be damaged by surgical instruments or bone cement particles. In order to guarantee the safe and reliable function, they have to remain intact. Therefore, these areas should be covered with surgical drapes during implantation.

5. Indications

Implantation of a prosthesis should be taken into consideration, if conservative measures or other surgical measures prove not to be appropriate. Even successfully implanted prostheses are inferior to the natural joint. Implant loosening, wear, and implant failure may necessitate revision surgery.

5.1 Indications REVISIO® RTM System

- Loss of the weight-bearing function of bones and adjacent joints
- Degenerative joint diseases (e.g. avascular necrosis, osteoarthritis, rheumatoid arthritis, requiring extended resection and reconstruction of the hip and/or knee joint).
- Significant bone loss and/or severe ligament instability.
- Correction of functional varus and valgus deformities
- Metastatic diseases (e.g. giant cell carcinoma, osteosarcoma, chondrosarcoma, bone tumours requiring extended resection and reconstruction of the proximal and/or diaphyseal femur and/or proximal tibia).
- Multiple revision arthroplasties of the knee and/or hip joint requiring reconstruction of the proximal, diaphyseal and/or distal femurs and/or proximal tibia.
- Traumatic injuries requiring extended replacement of the proximal, distal, diaphyseal or total femur and/or proximal tibia.
- In cases where alternative, less invasive methods do not seem to be promising

5.2 Specific indications REVISIO® RTM System

- Proximal femur replacement: large bone defects (tumours) in the proximal femur
- Diaphyseal replacement: Bone defects (tumours) in the diaphyseal femur
- Distal femur replacement: Bone defects (tumours) in the distal femur with partial replacement of the knee joint preserving the proximal tibia
- Distal femur replacement and proximal tibia replacement: Bone defects (tumours) in the distal femur and proximal tibia with total replacement of the knee joint.
- Total femur replacement: Bone defects (tumours) in the femur and tibial head with total

replacement of the knee joint.

5.3 Indications REVISIO® RTM Endofemoral Module

The intramedullary total femur is indicated in cases where load-stable reconstruction of the femur can be achieved with or without partial resection and existing tendon and muscle insertions as well as the patellar bearing can be preserved.

- Loss of the weight-bearing function of bones and adjacent joints
- Degenerative joint diseases (e.g. avascular necrosis, osteoarthritis, rheumatoid arthritis, requiring extended resection and reconstruction of the hip and/or knee joint).
- Significant bone loss and/or severe ligament instability.
- Correction of functional varus and valgus deformities
- Metastatic diseases (e.g. giant cell carcinoma, osteosarcoma, chondrosarcoma, bone tumours requiring extended resection and reconstruction of the proximal and/or diaphyseal femur and/or proximal tibia).
- Multiple revision arthroplasties of the knee and/or hip joint requiring reconstruction of the proximal, diaphyseal and/or distal femurs and/or proximal tibia.
- Traumatic injuries requiring extended replacement of the proximal, distal, diaphyseal or total femur and/or proximal tibia.
- Periprosthetic fractures
- Intraprosthetic fractures
- According to the indication for the REVISIO® RTM systems, whereas existing tendon and muscle insertions as well as the patellar bearing might be preserved.

Note:

Indication for cementless procedures requires sufficient bone quality. Especially in patients with cement incompatibility, cementless fixation provides a viable alternative. If cementless fixation is contraindicated, the prosthesis must be cemented.

6. Contraindications

6.1 General contraindications

- Acute or chronic, local or systemic, uncontrollable infection
- Insufficient bone quality or conditions jeopardising a stable fixation of the prosthesis
- Vascular, muscular, or neurological diseases compromising the affected area
- Any comorbidity that might affect bone regeneration or implant function (e.g. severe metabolic disorders)
- High physical strain (e.g. competitive sports, hard physical work)
- Incomplete skeletal growth

- Allergies to implant materials
- Morbid obesity or high weight gain
- Pregnancy and breastfeeding
- Lack of compliance
- Intake of high doses of cortisone or cytostatic drugs
- High tendency of falling or other circumstances that might lead to an increased fracture risk
- Cement intolerance in cemented prostheses

Contraindications comprise any circumstances making surgery inadvisable or impossible, such as coagulation disorders (also caused by indispensable anticoagulants), non-treatable high blood pressure, new myocardial infarction, cardio-vascular diseases, severe cardiac insufficiency, anaesthetic intolerance, or poor general condition of the patient.

Mental and neurological conditions present a contraindication in case of the patient's lack of willingness or ability to restrict movement.

There is no adequate research on the use during pregnancy and breastfeeding. Therefore, the implant should not be used in pregnant and nursing women.

6.2 Specific contraindications

a) REVISIO® RTM System

- Contraindications regarding metastatic diseases:
- Any concomitant condition that might lead to a local or non-local spread of tumours
- Pathological bone fractures
- Open infections
- Excessively fast progression of the disease

b) Endofemoral Replacement

- Acute purulent bone or soft tissue inflammations
- Cement ageing and specific antibiotics in the cement mixture

7. Potential adverse effects

The expected life expectancy of an implant can be impaired in young, very tall and obese patients and patients who are physically very active. Wrong selection or maladjustment of the

prosthesis, improper fixation, positioning or use despite of contraindication, can lead to early failure due to loosening, fracture, or wear.

Typical adverse effects of artificial joints may include:

- Implant loosening or displacement and wear
- Dislocation or loosening and fracture of components
- Infection
- Impingement
- Change of limb length
- Periarticular calcification
- Heterotopic ossification
- Hypersensitivity to alloy components
- Periprosthetic fractures
- Chronic pain
- Pulmonary embolism and venous thrombosis
- Nerve damage, wound healing disorder, and haematoma
- Cardiovascular disorders
- Limited joint function and mobility
- Limited joint loading and joint pain
- Changes of the bone tissue (atrophy, necrosis)

8. Patient information

Serial or batch number as well as item number of the implanted prosthesis should be documented in the patient's record and the implant pass the patient receives after surgery. Labels are included in the sterile packaged implants.

The patient must be informed about potential benefits and risks of the procedure and the expected effects resulting from individual physiological conditions, e.g. bone quality.

The patient must be informed about the proper handling of the implant and avoidance of harmful activities that might lead to implant loosening or fatigue as well as bone fractures. The patient should also be informed that the degree of motion prescribed during the informed consent discussion must be observed. Moreover, the patient should be informed that sports activities are possible to a limited extend and after consultation and individual risk assessment by the attending physician.

Patient information should also include the possible need for partial or total revision surgery, e.g. due to implant loosening as a consequence of mechanical wear or infection.

The patient must be informed that the body weight plays a decisive role in joint replacement and that heavy loads should be avoided.

However, the patient must be informed that regular physiotherapy should be carried out in accordance with the individual instructions of the attending physician.

The patient must be informed that diagnostic or follow-up magnetic resonance imaging (MRI) can cause interactions with the implant.

The patient should be informed that the use of cups or asymmetrical cup inlays can result in a limited range of motion (ROM).

Finally, the patient should be advised in which circumstances the attending physician should be contacted.

9. Implant materials

Depending on the individual indication, joint replacement components are implanted with or without bone cement (PMMA). The following implant materials are used:

- Cobalt alloy (CoCrMo) manufactured according to ISO 5832-4 or ISO 5832-12, respectively
- UHMWPE – Ultra-high molecular weight polyethylene according to ISO 5834-1, -2
- Calcium phosphate (CaP) coatings
- Titanium-Niobium [TiNb] or TiNb-ON coating of CoCrMo components for patients suffering from allergic reactions to nickel, cobalt, chromium, molybdenum
- Titanium alloy (TiAl6V4), produced according to ISO 5832-3

Note:

Please follow the instructions regarding the mechanical compatibility of the implant materials and consider potential allergic reactions to one or multiple components of these implant materials (see Contraindications).

10. Clinical follow-up

Patients with an artificial joint require continuous follow-up by the surgeon or another competent physician. In the event of an unforeseen early implant loosening, material weakening or similar signs of potential complications, appropriate measures can be taken in time. If a revision procedure should become necessary, the chances for success are significantly higher in case of early identification, especially if the patient's bone quality at the implant-bone interface has not yet

deteriorated significantly.

The patient should be advised to inform the surgeon or attending physician about any changes relating to the artificial joint.

Regular clinical follow-up examinations are recommended. If this is not possible, the surgeon or attending physician should be provided with a control X-ray for assessment at least once a year.

AQ Implants GmbH provides appropriate instruments for revision procedures. Please refer to the specific instructions in our user information and surgical techniques. Upon request, our product specialists will provide support during revision surgery at your convenience.

Caution:

In case of any queries relating revision procedures, please contact the implant manufacturer.

Caution:

Exposure to strong magnetic or electromagnetic fields commonly used in magnetic resonance imaging (MRI) can have adverse effects on the implant behaviour. For the latest safety information, please contact the implant manufacturer.

11. Compatibility

We guarantee compatibility of our products only for our CE-marked products and products we have approved for combination. Any combination of implants of AQ Implants GmbH with components of other manufacturer is excluded for reasons of product safety and liability.

The modular Revision and tumour system REVISIO® RTM can be combined with metal heads and with BIOLOX® Delta and BIOLOX® Forte ceramic heads by Ceramtec. In case of any queries please contact us. Please have the item-no. and ident-no. (lot-no., where appropriate) of the respective products ready.

The REVISIO® RTM Cone Adapter provides the basis for a universal and versatile implant system for the treatment of complex defects. Only in specific cases, the REVISIO® RTM system is indicated for adolescents or children. It is important to note that the patient's bone conditions can change significantly during growth which might compromise stable primary fixation of the implant. Given this situation, supportive measures should be taken into consideration. In any case, follow-up x-rays should be taken at regular tight intervals to identify any changes that might pose a risk to the implant and take appropriate measures immediately.

Please ensure stable fixation of the adapter stem over the entire implant length in a closed cement mantle, particularly in small cemented stems, in order to ensure maximum fatigue strength for natural loading.

Depending on the complexity of treatment (bridging of large defects), the biomechanical loading of the adapter stems can increase significantly, e.g. due to the patient's high body weight or increased activity level. In those circumstances, in particular in small stem sizes, appropriate measures should be taken, e.g. partial weight-bearing.

Any non-observance of these recommendations can lead to fracture of the implant in extreme cases, in particular in small stem sizes.

User information for REVISIO® RTM Cone Adapter (adapter hip stem GHE, GHEs and REVISIO® RTM Trochanteric Module Endofemoral)

The modular design of hip prostheses with REVISIO® RTM Cone Adapter allows for patient-individual reconstruction of the hip joint. Due to an increased lever arm, loading of the hip prostheses is significantly increased in REVISIO® RTM Cone Adapters with an angle of >0° and lateralisation. This effect is enhanced with increasing adapter angle and small stem sizes. For the REVISIO® RTM Trochanteric Module, REVISIO® RTM Cone Adapter with 0, 5, 10° and lateralised (size 4) can be used.

The REVISIO® RTM cone adapter may only be combined with femoral heads (28, 32, 36mm) in lengths S/ML.

Any combination of REVISIO® RTM cone adapter with REVISIO® RTM Extensio modules is not permitted!

In case of any deviating combination due to patient-individual factors and thus clinically required indication, utmost care and diligence are recommended, including detailed patient information and observance of all required safety precautions (e.g. partial weight-bearing).

Especially in small endoprostheses or small stem sizes of the REVISIO® RTM Trochanteric Modul, respectively, as well as obese and very active patients, the use of walking aids should be recommended. Any application adapter angles of >10° is only suitable upon request and as a stopgap in extreme situations in compliance with the safety precautions.

Please refer to the specific instructions for cone adapters in our surgical technique for these hip prosthesis systems or contact our Customer Service.

User information for the use of extension adapters

If extension adapters are used, an additional neck lengthening can be achieved. Thus, they should be used in patients with hip endoprosthesis with dislocation tendencies that cannot be

prevented with extra large balls. Extension adapters increase the load of the hip stem significantly, as they extend the lever arm of the point of load application.

Therefore, the following points should be observed:

- The extension adapter should not be used with small stems (size 1) to avoid overloading the component.
- Adapters of sizes 2 to 4 should not be used in obese patients to avoid implant fractures.
- The extension adapter should not be combined with the REVISIO® RTM cone adapter and balls with a large neck length (XL-XXXL).

In case of any deviating combination due to patient-individual factors and thus clinically required indication, utmost care and diligence are recommended, including detailed patient information and observance of all required safety precautions.

12. Instrument care – special precautions

Proper handling, decontamination, including pre-rinsing, washing, rinsing and sterilisation, storage, and use of instruments are decisive for a long and adequate durability of all surgical instruments. However, even with correct use, maintenance and care, life expectancy of the instruments is limited. This applies, in particular, to cutting instruments and force-controlled instruments. These instruments and their components are often exposed to high loads and/or impaction forces. Defects can occur under such conditions, especially if the instrument is corroded, damaged, notched, or scratched.

Please replace any instruments that have become unusable. Examine all instruments carefully before using them. Never use an instrument that is damaged or worn, or has an illegible laser engraving or a cutting instrument with blunt cutting edges.

Note:

Please follow the general instructions on the maintenance and care of surgical instruments. Please refer to the „Reprocessing guidelines for the reprocessing of reusable surgical instruments“ of AQ Implants GmbH for detailed information. In case of any doubt, please contact us.

13. Symbols

Appendix 2 contains the symbols used on the package labels.

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AQ Implants GmbH tarafından kullanılan semboller, 2 ekinde görülebilir.