**MDD regulation**

**Gebrauchsanweisung REVISIO® RTM**  
**Instructions for use REVISIO® RTM**

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## 1. Instructions for Use Revision and Tumour System

The present Instructions for Use applies to the REVISIO® RTM revision and tumour system. Appendix 1 contains an implant overview including all relevant item numbers.

Improper implantation of the prosthesis can lead to failure, implant failure, patient injury, bone loss, or unsatisfactory aesthetic results. The products may only be used by trained physicians and surgeons who are familiar with the implant system. Use of these products requires specialised knowledge and skills in endoprosthesis/surgical implantology. AQ Solutions GmbH offers specific training and technical consultation on the use of its products and provides detailed information for appropriate implant selection as well as preoperative planning and surgical technique. Our products are sold only to physicians and healthcare professionals or on their behalf.

Technical consultation on our products will be performed verbally or in writing, by electronic media, and/or by demonstration. Consultation reflects the latest state of science and technology at the time of the product launch. It does not relieve the user from the obligation to personally verify suitability of the product in terms of intended purpose, indications, and procedures. Handling and implantation of the product are beyond our control and are the responsibility of the user.

## 2. Product description and intended use

The REVISIO® RTM system was developed for complex applications in tumour surgery. Due to the high modularity of the REVISIO® RTM system, the individual components of the proximal and distal femur as well as the diaphysis up to the complete femur can be replaced. The REVISIO® RTM system also features an arthrodesis spacer for joint arthrodesis of the knee joint as well as the REVISIO® RTM condylar replacement to maintain the mobility of the knee joint. The REVISIO® RTM endofemoral module which is introduced through the femur connecting hip and knee implant represents an expansion of the REVISIO® RTM system. It was developed to minimise bone loss in the femur by preserving the diaphyseal femur.

The modular implant components are available in various diameters and lengths ensuring an adequate reconstruction of the individual defect.

The individual modules are connected by means of tapered clamping and additional fixation via a specific screw system.

The combination of implant components depends on the indication and is the responsibility of the surgeon.

### 3. Handling

Keep the implant in the unopened original packaging until it is used. When unpacking the implant, please ensure that the protective packaging, i.e. cardboard and peel pouch is unopened and intact. Any damage may impair the sterility of the implant. In the event of a broken seal or damaged cardboard and peel pouch, it must be assumed that the contents are no longer sterile.

**Caution:**

If the seal is broken, the product cannot be returned.

When unpacking the implants, please verify the details on the packaging against the package contents (item number and size) as well as expiry date. Do not implant components after the expiry date stated on the package.

**Caution:**

Make sure that all aseptic standards are observed during removal of sterile implant components.

Take appropriate measures to avoid that the implant comes into contact with objects or substances that might damage its surface. This applies particularly for coated implants (Titanium-Niobium [TiNb] etc.) e.g. for the contact with body fluids (sweat, blood, synovia, etc.). For ceramic implants this also applies for the direct contact with skin. Please note the enclosed user information and surgical techniques. Visually inspect implant and coating (where applicable) for damage prior to implantation. Damaged implants or implant components must not be implanted. Protective caps or similar devices for the protection of sensitive implant areas must be removed immediately before use. Implants must not be processed mechanically, altered or modified, unless this is explicitly required by the surgical technique or the manufacturer has given his written consent.

**Caution:**

Products that have been implanted must not be reused. The same applies for implants that came into contact with foreign body fluids.

Removal, handling and return of revised implants for scientific evaluation must be performed in

accordance with ISO 12891-1. Please observe the general information of ISO 8828 on the care and handling of orthopaedic implants! For specific items, please observe the manufacturer's instructions (handling, sterilisation, compatibility, etc.).

#### 4. Packaging and sterility

All implants are packed individually in sterile protective packages. Sterile implants have been sterilised with a minimum gamma radiation dose of 25 kGy.

 **Caution:** Any implants labelled „sterile“ are only sterile, if the package has not been opened.

 **Caution:** Implants must not be resterilised. Any reprocessing of non-implanted components with opened packages or expired sterile date must only be performed by the manufacturer, as several validated processes have to be performed again.

Implants intended for use with bone cement: C (cemented)

#### 5. Surgical planning

Preoperative planning is based on the radiographic evaluation. Appropriate X-ray templates are available. Moreover, all products of AQ Solutions GmbH are available in the digital planning system mediCAD®.

For the preparation of the host bone and adjustment, repositioning as well as implantation of the implants, specific instruments and trial implants corresponding to the dimensions of the implants are available. Please don't use other trial instruments or implants. The appropriate surgical technique is selected in accordance with the latest state of science and technology. In addition, AQ Solutions GmbH provides surgical techniques upon request which are constantly revised in view of new findings and understanding.

In order to achieve optimal primary fixation as a prerequisite for long-term stability of the prosthesis, the following points should be observed:

- Preservation of adequate bone structure / wall thickness during the preparation of the host bone
- Selection of the appropriate implant size
- Adequate cleaning of the host bone prior to implantation
- Application of moderate force during implantation to avoid bone fissures or implant damage.
- Compatibility of the individual components
- Observance of surgical technique, user notes, and warnings
- Clean and dry surfaces of tapered connections

- Firm coupling of tapers achieved by using insertion and impaction instruments
- Avoid damage to the components
- Avoid foreign particles, such as bone or tissue debris or cement particles between functional surfaces (e.g. tapers, holes, articular surfaces)
- Use of locking screws in screw-retained components (e.g. fixation of the iliac peg)

The tapers of the prostheses can be damaged by surgical instruments or bone cement particles. In order to guarantee the safe and reliable function, they have to remain intact. Therefore, these areas should be covered with the supplied taper covers and surgical drapes during implantation. Please ensure stable fixation of the RTM stem over the entire implant length in a closed cement mantle, particularly in small cemented stems, in order to ensure maximum fatigue strength for natural loading.

Depending on the complexity of treatment (bridging of large defects), the biomechanical loading of the adapter stems can increase significantly, e.g. due to the patient's high body weight or increased activity level. In those circumstances, in particular in small stem sizes, appropriate measures should be taken, e.g. partial weight-bearing.

Any non-observance of these recommendations can lead to fracture of the implant in extreme cases, in particular in small stem sizes.

#### **Cautions:**

- Application of moderate force during implantation to avoid bone fissures or implant damage
- Avoid foreign particles, such as bone or tissue debris or cement particles between functional surfaces (e.g. tapers, articular surfaces)
- The following notes must be observed for the process for implanting cemented implant components:
  - Preparation and application of bone cement according to instructions provided by the respective manufacturer
  - Sufficient cleaning and drying of the bone site prior to implantation
  - Positioning of a medullary canal stopper in the marrow and insertion of the bone cement or alternatively application of the cement directly onto the implant
  - Prompt and complete removal of leftover remains of cement
- Insufficient securing of the cone can lead to screw loosening or to loosening of the cone, which could require a further corrective operation.
- It must be ensured that the spacer components with the ventral bone are positioned along the tibia, to allow the coupling of the spacer using the locking screw.

- We guarantee the compatibility of our products only in connection with our own products and the products from other manufacturers we have approved for the combination. The combination of AQ Solutions implants with not approved products from other manufacturers is excluded for reasons of product safety and product liability. Please refer to the note on the compatibility of products from AQ Solutions GmbH in our OP-manual.
- The metal surfaces must be free of any residues. In some cases the components may need to be cleaned and dried several times! Absolutely all surfaces must be entirely clean.
- For later revisions, it is important to note that the threaded opening of the cone adaptor faces the head.
- When using a ceramic head, please note the manufacturers specific instructions for implantation and fixation.
- If the defect size requires the use of an Extensio module, please note that, due to a minimum component length, the next resection level must be 100 mm. The reconstruction can then be carried out in 10 mm steps.
- An inadequate securing of the cone can lead to screw loosening or to loosening of the cone, which could necessitate a further revision operation.
- The inner and outer cones must be free of contaminant and liquid.
- AQ Solutions GmbH supplies cement-based and cement-less acetabulum and alternatively, bipolar heads or large heads. The REVISIO® RTM system is compatible with all AQ Solutions GmbH implants for acetabulum treatment. The combination of REVISIO® RTM System and REVISIO® M Cup System is particularly suitable for restorations with large defects.
- The cement mantle should fully enclose the implant in the fixation area. An incomplete cement mantle can lead to early loosening of the implant.
- Bearing surfaces and tapers of the prostheses can be damaged by surgical instruments or bone cement particles. In order to guarantee the safe and reliable function, they have to remain intact. Therefore, these areas should be covered with surgical drapes during implantation

### **Caution for condyle replacement**

- Please note that the fixing screw is placed loosely in the metal support and must be removed before implantation.

## **6. Indications**

### **6.1 Indication REVISIO® RTM System**

REVISIO® RTM-System is indicated for treatment of lower extremities where extensive bone and/or joint reconstruction is required in proximal and/or diaphyseal and/or distal femur and/or pro-

ximal tibia in revision or tumour cases.

Implantation of a prosthesis should be taken into consideration, if conservative measures or other surgical measures prove not to be appropriate. The REVISIO® RTM-System is indicated, in cases where alternative, less invasive methods do not seem to be promising.

The indications can result from any of the following medical conditions:

- Primary and secondary bone tumors resulting in large bone defects
  - in the proximal or diaphyseal region of the femur
  - in the distal region of the femur with partial replacement of the knee joint while preserving the proximal tibia
  - in the distal region of the femur and proximal tibia with complete replacement of the knee joint
  - throughout the femur and tibial head with total replacement of the knee joint
  - including metastatic diseases (e.g. giant cell carcinoma, osteosarcoma, chondrosarcoma, bone tumours requiring extended resection and reconstruction of the proximal and/or diaphyseal femur and/or proximal tibia)
- Loss of the weight-bearing function of bones and adjacent joints
- Degenerative joint diseases (e.g. avascular necrosis, osteoarthritis, rheumatoid arthritis, requiring extended resection and reconstruction of the hip and/or knee joint)
- Significant bone loss and/or severe ligament instability
- Multiple revision arthroplasties of the knee and/or hip joint requiring reconstruction of the proximal, diaphyseal and/or distal femurs and/or proximal tibia
- Traumatic injuries requiring extended replacement of the proximal, distal, diaphyseal or total femur and/or proximal tibia
- Periprosthetic fractures
- Intraprosthetic fractures

## 6.2 Indications REVISIO® RTM femoral metal heads

The metal heads are indicated in conjunction with the REVISIO® RTM system for proximal femoral replacement, total endofemoral replacement or total femoral replacement.

## 7. Contraindications

- Acute or chronic, local or systemic, unmanageable infection
- Deficient bone quality or bone conditions that endanger a stable anchorage of the prosthesis
- Vascular, muscle and nerve diseases that endanger the affected area

- Any accompanying disease that can impair bone regeneration or the implant's function (e.g. the most severe metabolic disorders)
- High physical strain (e.g. competitive sports, heavy physical labour)
- Incomplete skeletal growth
- Allergies to components of implant materials
- Extreme obesity or severe weight gain
- Pregnancy and breastfeeding
- Lack of compliance
- Taking high doses of cortisone or cytostatic drugs
- Epilepsy or other circumstances which could lead to accidents with an increased risk of fractures
- Cement intolerances with cemented implants
- Soft tissue inflammations
- Contraindications regarding metastatic diseases:
  - Any concomitant condition that might lead to a local or non-local spread of tumors
  - Pathological bone fractures
  - Open infections
  - Excessively fast progression of the disease

Contraindications also include those that generally exclude surgery, such as coagulation disorders (also due to anticoagulants that are indispensable), untreatable high blood pressure, recent heart attack, severe cardiovascular diseases, severe cardiac insufficiency, intolerance to anesthesia or poor general condition of the patient. Mental and neurological conditions present a contraindication in case of the patient's lack of willingness or ability to restrict movement. These contraindications are not related to the implant itself and are at the sole discretion of the physician in charge based, among other aspects, on the presentation of the patient, the specific diagnosis/es, diagnostic findings, pathophysiologic reasoning, and personal experience.

## 8. Potential adverse effects

The expected life expectancy of an implant can be impaired in young, very tall and obese patients and patients who are physically very active. Wrong selection or maladjustment of the prosthesis, improper fixation, positioning or use despite of contraindication, can lead to early failure due to loosening, fracture, or wear.

Typical adverse effects of artificial joints may include:

- Contamination of bone marrow space with bone cement if bone marrow space stop is not used
- Limited implant survival if patient is severely obese
- Aseptic implant loosening
- Chronic pain due to inflammation
- Periprosthetic infection and septic loosening
- Limited joint function and mobility due high inherent weight
- Hypersensitivity to alloy components
- Luxation
- Periprosthetic fractures
- Impingement
- Dislocation, loosening or fracture of components
- Change of limb length
- Changes of the bone tissue (atrophy, necrosis)
- Nerve damage, wound healing impairment, and haematoma
- Failure of the hinge mechanism

## 9. Patient information

Serial or batch number as well as item number of the implanted prosthesis should be documented in the patient's record and the implant pass the patient receives after surgery. Labels are included in the sterile packaged implants.

The patient must be informed about potential benefits and risks of the procedure and the expected effects resulting from individual physiological conditions, e.g. bone quality.

The patient must be informed about the proper handling of the implant and avoidance of harmful activities that might lead to implant loosening or fatigue as well as bone fractures. The patient should also be informed that the degree of motion prescribed during the informed consent discussion must be observed. Moreover, the patient should be informed that sports activities are possible to a limited extend and after consultation and individual risk assessment by the attending physician.

Patient information should also include the possible need for partial or total revision surgery, e.g. due to implant loosening as a consequence of mechanical wear or infection.

The patient must be informed that the body weight plays a decisive role in joint replacement and that heavy loads should be avoided.

However, the patient must be informed that regular physiotherapy should be carried out in accordance with the individual instructions of the attending physician.  
The patient must be informed that diagnostic or follow-up magnetic resonance imaging (MRI) can cause interactions with the implant.

The patient should be informed that the use of cups or asymmetrical cup inlays can result in a limited range of motion (ROM).  
Finally, the patient should be advised in which circumstances the attending physician should be contacted.

Patient information should also include the possible need for partial or total revision surgery, e.g. due to implant loosening as a consequence of mechanical wear or infection.

The patient must be informed that the body weight plays a decisive role in joint replacement and that heavy loads should be avoided.

However, the patient must be informed that regular physiotherapy should be carried out in accordance with the individual instructions of the attending physician.

The patient must be informed that diagnostic or follow-up magnetic resonance imaging (MRI) can cause interactions with the implant.

The patient should be informed that the use of cups or asymmetrical cup inlays can result in a limited range of motion (ROM).

Finally, the patient should be advised in which circumstances the attending physician should be contacted.

## 10. Implant materials

Depending on the individual indication, joint replacement components are implanted with bone cement (PMMA). The following implant materials are used:

- Cobalt alloy (CoCrMo) manufactured according to ISO 5832-4 or ISO 5832-12, respectively
- UHMWPE – Ultra-high molecular weight polyethylene according to ISO 5834-1, -2
- Titanium-Niobium [TiNb] or TiNb-ON coating of CoCrMo components for patients suffering from allergic reactions to nickel, cobalt, chromium, molybdenum

- Titanium alloy (TiAl6V4), produced according to ISO 5832-3

**Note:**

Please follow the instructions regarding the mechanical compatibility of the implant materials and consider potential allergic reactions to one or multiple components of these implant materials (see Contraindications).

## 11. Clinical follow-up

Patients with an artificial joint require continuous follow-up by the surgeon or another competent physician. In the event of an unforeseen early implant loosening, material weakening or similar signs of potential complications, appropriate measures can be taken in time. If a revision procedure should become necessary, the chances for success are significantly higher in case of early identification, especially if the patient's bone quality at the implant-bone interface has not yet deteriorated significantly.

The patient should be advised to inform the surgeon or attending physician about any changes relating to the artificial joint.

Regular clinical follow-up examinations are recommended. If this is not possible, the surgeon or attending physician should be provided with a control X-ray for assessment at least once a year.

AQ solutions GmbH provides appropriate instruments for revision procedures. Please refer to the specific instructions in our user information and surgical techniques. Upon request, our product specialists will provide support during revision surgery at your convenience.

**Caution:**

In case of any queries relating revision procedures, please contact the implant manufacturer.

**Caution:**

Exposure to strong magnetic or electromagnetic fields commonly used in magnetic resonance imaging (MRI) can have adverse effects on the implant behavior. For the latest safety information, please contact the implant manufacturer.

## 12. Compatibility

We guarantee compatibility of our products only for our CE-marked products and products we have approved for combination. Any combination of implants of AQ Solutions GmbH with components of other manufacturer is excluded for reasons of product safety and liability.

The modular Revision and tumour system REVISIO® RTM can be combined with metal heads and with BIOLOX Delta and BIOLOX Option ceramic heads. In case of any queries please contact us. Please have the item-no. and ident-no. (lot-no., where appropriate) of the respective products ready.

### **User information for REVISIO® RTM Cone Adapter and REVISIO® RTM Trochanteric replacement**

The modular design of hip prostheses with REVISIO® RTM Cone Adapter allows for patient-individual reconstruction of the hip joint. Due to an increased lever arm, loading of the hip prostheses is significantly increased in REVISIO® RTM Cone Adapters with an angle of  $>0^\circ$  and lateralisation. This effect is enhanced with increasing adapter angle and small stem sizes. For the REVISIO® RTM Trochanteric replacement, REVISIO® RTM Cone Adapter with 0, 5, 10° and lateralised (size 4) can be used.

The REVISIO® RTM cone adapter may only be combined with femoral heads (28, 32, 36mm) in lengths S/ML.

Any combination of REVISIO® RTM cone adapter with REVISIO® RTM Extensio modules is not permitted!

In case of any deviating combination due to patient-individual factors and thus clinically required indication, utmost care and diligence are recommended, including detailed patient information and observance of all required safety precautions (e.g. partial weight-bearing).

Especially in small endoprotheses or small stem sizes of the REVISIO® RTM Trochanteric replacement, respectively, as well as obese and very active patients, the use of walking aids should be recommended. Any application adapter angles of  $>10^\circ$  is only suitable upon request and as a stopgap in extreme situations in compliance with the safety precautions.

Please refer to the specific instructions for cone adapters in our surgical technique for these hip prosthesis systems or contact our Customer Service.

**User information for the use of extension adapters**

If extension adapters are used, an additional neck lengthening can be achieved. Thus, they should be used in patients with hip endoprosthesis with dislocation tendencies that cannot be prevented with extra large balls. Extension adapters increase the load of the hip stem significantly, as they extend the lever arm of the point of load application.

- The extension adapter should not be used with small stems (size 1) to avoid overloading the component.
- Adapters of sizes 2 to 4 should not be used in obese patients to avoid implant fractures.
- The extension adapter should not be combined with the REVISIO® RTM cone adapter.

In case of any deviating combination due to patient-individual factors and thus clinically required indication, utmost care and diligence are recommended, including detailed patient information and observance of all required safety precautions.

**13. Symbols**

Appendix 2 contains the symbols used on the package labels.

Bezeichnung • Designation	Artikelnummer • Item number		
Metallkopf Metal head	16000128 16000228 16000328	16000132 16000232 16000332	16000136 16000236 16000336
REVISIO® RTM Konusadapter REVISIO® RTM Cone adapter	60046011 - 60046017 60046019 - 60046028		
REVISIO® RTM Trochantermodul REVISIO® RTM Trochanteric replacement	60050007 - 60050010 60045003 - 60045007 60045103 - 60045107 60045203 - 60045207 60045303 - 60045307		
REVISIO® RTM Extensiomodul REVISIO® RTM Extensio Module	60055030 60055060 60055160 60056050 60056120	60055040 60055080 60056030 60056060 60056160	60055050 60055120 60056040 60056080
REVISIO® RTM Connector REVISIO® RTM Connector	60057001 - 60057014		
REVISIO® RTM Sicherungsschraube REVISIO® RTM Locking screw	60058001 - 60058004 60058008 60058016	60058010 60058020	60058012
REVISIO® RTM Endofemoral-Modul REVISIO® RTM Endofemoral Module	60059100 - 60059125 60059200 - 60059225		

Bezeichnung • Designation	Artikelnummer • Item number
REVISIO® RTM Stiel C REVISIO® RTM Stem C	60048011 - 60048016 60048111 - 60048116 60048211 - 60048216 60048311 - 60048316 60048411 - 60048416 60048511 - 60048516
REVISIO® RTM Femoral-Tibial-Stiel REVISIO® RTM Femoral tibial stem	60067001 - 60067010 60062201 - 60062206 60062301 - 60062306
REVISIO® RTM Arthrodesen-Spacer REVISIO® RTM Arthrodesis spacer	60063000    60063025    60063050 60063125    60063150    60063300 60063325    60063350    60063425 60063450
PE-PATELLA ALL-POLY-PATELLA	27700010    27700020    27700030
REVISIO® RTM Kondylenersatz REVISIO® RTM condylar replacement	60049001 - 60049027
REVISIO® RTM Tibiakopf REVISIO® RTM Tibial head	60054001 - 60054003 60054101 - 60054103
REVISIO® RTM Femurkasten REVISIO® RTM Femoral box	60060001 - 60060004 60072101 - 60072106 60072001 - 60072006 60060101 - 60060104
REVISIO® RTM Femur-Gelenkteil REVISIO® RTM Femoral joint module	60061001 - 60061004 60061101 - 60061104

Bezeichnung • Designation	Artikelnummer • Item number		
REVISIO® RTM Tibia-Gelenkteil REVISIO® RTM Tibial component	60064001 - 60064006		
REVISIO® RTM Tibia-Plateau REVISIO® RTM Tibial Plate	60065101 - 60065106 60074101 - 60074105 60074001 - 60074005		
REVISIO® RTM Tibia-Augmentat REVISIO® RTM Tibial augment	60066008 60066016 60066110 60066120 60066212 60066308 60066316	60066010 60066020 60066112 60066208 60066216 60066310 60066320	60066012 60066108 60066116 60066210 60066220 60066312
REVISIO® RTM Femur-Gelenkteil REVISIO® RTM Femoral component	60071001 - 60071006 60071101 - 60071106		
REVISIO® RTM Tibia-Einsatz REVISIO® RTM Tibia Insert	60073108 60073114 60073212 60073310 60073410 60073414	60073110 60073208 60073214 60073312 60073412	60073112 60073210 60073308 60073408 60073314
REVISIO® RTM Tibiastiel REVISIO® RTM Tibial stem	60075001 - 60075005 60075101 - 60075105		

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