MDD regulation

Gebrauchsanweisung Acetabulum
Instructions for use Acetabulum

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Instructions for use Acetabulum

The present Instructions for Use Acetabulum apply to the following implants:

- **REVISIO® S Acetabular Cup**
- **REVISIO® M Acetabular Cup Standard**
- **REVISIO® M Acetabular Cup with flange**
- **REVISIO® PE Inlays Standard**
- **REVISIO® PE Inlays Asymmetric 10°/20° and Constrained**
- **Accessories**

Appendix 1 contains an implant overview including all relevant item numbers.

Improper implantation of the prosthesis can lead to failure, implant failure, patient injury, bone loss, or unsatisfactory aesthetic results. The products may only be used by trained physicians and surgeons who are familiar with the implant system. Use of these products requires specialised knowledge and skills in endoprosthetic/surgical implantology. AQ Solutions GmbH offers specific training and technical consultation on the use of its products and provides detailed information for appropriate implant selection as well as preoperative planning and surgical technique. Our products are sold only to physicians and healthcare professionals or on their behalf.

Technical consultation on our products will be performed verbally or in writing, by electronic media, and/or by demonstration. Consultation reflects the latest state of science and technology at the time of the product launch. It does not relieve the user from the obligation to personally verify suitability of the product in terms of intended purpose, indications, and procedures. Handling and implantation of the product are beyond our control and are the responsibility of the user.

1. Product description and intended use

REVISIO® S and M Acetabular Cup are modular acetabular systems for cementless fixation in the acetabulum. REVISIO® S and M Acetabular Cup have a 3D macroporous surface and are anchored press-fit in the acetabulum. They are combined with specific PE-Inlays that have been developed for these systems.

The combination of implant components depends on the indication and is the responsibility of the surgeon.

2. Handling

Keep the implant in the unopened original packaging until it is used. When unpacking the implant, please ensure that the protective packaging, i.e. cardboard and peel pouch, is unopened and intact. Any damage may impair the sterility of the implant. In the event of a broken seal or damaged cardboard and peel pouch, it must be assumed that the contents are no longer sterile.

Caution:

If the seal is broken, the product cannot be returned.

When unpacking the implants, please verify the details on the packaging against the package contents (item number and size) as well as expiry date. Do not implant components after the expiry date stated on the package.

Caution:

Make sure that all aseptic standards are observed during removal of sterile implant components.

Take appropriate measures to avoid that the implant comes into contact with objects or substances that might damage its surface. This applies particularly for coated implants (Titanium-Niobium [TiNb] etc.) e.g. for the contact with body fluids (sweat, blood, synovia, etc.). Please note the enclosed user information and surgical techniques. Visually inspect implant and coating (where applicable) for damage prior to implantation. Damaged implants or implant components must not be implanted. Protective caps or similar devices for the protection of sensitive implant areas are removed immediately before use. Implants must not be processed mechanically, altered or modified, unless this is explicitly required by the surgical technique or the manufacturer has given his written consent.

Caution:

Products that have been implanted must not be reused. The same applies for implants that came into contact with foreign body fluids.

Removal, handling and return of revised implants for scientific evaluation must be performed in accordance with ISO 12891-1. Please observe the general information of ISO 8828 on the care and handling of orthopaedic implants! For specific items, please observe the manufacturer's instructions (handling, sterilisation, compatibility, etc.).

3. Packaging and sterility

All implants are packed individually in sterile protective packages. Sterile implants have been sterilised with a minimum gamma radiation dose of 25 kGy.

Caution: Any implants labelled „sterile“ are only sterile, if the package has not been opened.

Caution: Implants must not be resterilised. Any reprocessing of non-implanted components with opened packages or expired sterile date must only be performed by the manufacturer, as several validated processes have to be performed again.

4. Surgical planning

Preoperative planning is based on the radiographic evaluation. Appropriate X-ray templates are available. Moreover, all products of AQ Solutions GmbH are available in the digital planning system medicAD®.

For the preparation of the host bone and adjustment, repositioning as well as implantation of the implants, specific instruments and trial implants corresponding to the dimensions of the implants are available. Please don't use other trial instruments or implants.

The appropriate surgical technique is selected in accordance with the latest state of science and technology. In addition, AQ Solutions GmbH provides surgical techniques upon request which are constantly revised in view of new findings and understanding.

In order to achieve optimal primary fixation as a prerequisite for long-term stability of the prosthesis, the following points should be observed:

- Preservation of adequate bone structure / wall thickness during the preparation of the host bone
- Selection of the appropriate implant size
- Adequate cleaning of the host bone prior to implantation
- Application of moderate force during implantation to avoid bone fissures or implant damage.
- Compatibility of the individual components
- Observance of surgical technique, user notes, and warnings
- Avoid damage to the components
- Use of locking screws in screw-retained components (e.g. fixation of the iliac peg)

Caution:

Bearing surfaces and tapers of the prostheses can be damaged by surgical instruments. In order to guarantee the safe and reliable function, they have to remain intact. Therefore, these areas should be covered with surgical drapes during implantation.

⚠ Caution:

Indication for cementless procedures requires sufficient bone quality. Especially in patients with cement incompatibility, cementless fixation provides a viable alternative.

⚠ Caution:

Application of moderate force during implantation to avoid bone fissures or implant damage

⚠ Caution:

Avoid foreign particles, such as bone or tissue debris between functional surfaces (e.g. tapers, articular surfaces)

⚠ Caution:

Excess osteophytes must be completely removed.

⚠ Caution:

The peg's length must not exceed the drilled hole's depth!

⚠ Caution:

To avoid damage to soft tissue, the screw length must be selected so that the screws do not extend beyond the drilling depth!

⚠ Caution:

The screw head must be completely sunk into the screw hole so as not to interfere with the clamping of the inlay!

The primary application of asymmetric inlays is recommended in exceptional cases, where there is an associated risk of dislocation, but re-positioning the cup is not feasible for stability reasons. The elevation of the PE Inlay must be positioned in such a way that it helps counteract the dislocation tendency. Prior to fixation it must be ensured that the elevated rim is in its predefined position.

The use of asymmetric inlays or acetabular cups, respectively, can result in a significant reduction of the range of motion (ROM). Especially the use of a 20° asymmetric inlay should thus be carefully considered.

**Caution:**

Use of PE Inlays with dislocation protection will restrict the postoperative mobility of the hip joint.

Caution: For PE Inlay Asymmetrical, the position of the superelevation is generally freely selectable. However, it is recommended to select the position of the superelevation in such a way that increased luxation security is achieved without femoro-acetabular impingement. Usually the superelevation is positioned in the area superior to dorsal acetabulum rim. With superior positioning the cup inclination angle is reduced, with dorsal positioning the cup inclination angle is increased.

Caution: Constrained System: Large amplitudes of movement will create a leverage effect on the locking mechanism which might damage the locking mechanism.

5. Indications

5.1 REVISIO® S Acetabular Cup

- Osteoarthritis of the hip of various etiologies
- Rheumatoid arthritis of the hip
- Pseudo-arthritis following femoral fractures
- Post-traumatic osteoarthritis of the hip
- Dysplastic hip joint
- Avascular necrosis of the femoral head
- Protrusio acetabuli
- Acetabular revision with smaller cavitory defects (Paprosky type I)
- The REVISIO® S hip cup screw fixation is specifically indicated in any of the above conditions when smaller cavitory defects (Paprosky type I) are present.

5.2 REVISIO® M Acetabular Cup Standard

Indication for an acetabular implant comprising the following disease entities associated with acetabular bone defects classified as Paprosky type IIa and IIb defect:

- Revision with acetabular bone defect situation
- Structural damage to the hip joint
- Tumour resections in the acetabular area
- Congenital hip dysplasia
- Osteoarthritis of the hip of various etiology

5.3 REVISIO® M Acetabular Cup with flange

Indication for an acetabular implant comprising the following disease entities associated with acetabular bone defects:

- Acetabular defects of type IIc to type IIIb according to Paprosky
- Revisions with acetabular bone defects of type IIc to type IIIb according to Paprosky
- Structural damage of the hip joint with resulting bone defects of type IIc to type IIIb according to Paprosky
- Osteoarthritis of the hip of various etiology with acetabular defects of type IIc to type IIIb according to Paprosky
- Tumour resections in the acetabular area resulting in bone defects of type IIc to type IIIb according to Paprosky

REVISIO® M acetabular cup with flange:

Paprosky type IIc and type IIIa defect classification

REVISIO® M acetabular cup with flange and fixation peg:

Paprosky type IIIa and type IIIb defect classification

5.4 REVISIO® PE-Inlays Standard

The REVISIO® PE-Inlays Standard are indicated in conjunction with the REVISIO® acetabular cups in patients with a normal risk of postsurgical dislocation

5.5 REVISIO® PE-Inlays Asymmetric 10°/20° and Constrained

The REVISIO® PE-Inlay 10° or 20° asymmetric, and the PE-Inlay constrained are indicated in conjunction with the REVISIO® acetabular cups with the following specific indications:

REVISIO® PE-Inlay 10° or 20° Asymmetric

- Postsurgical dislocation of the hip joint endoprosthesis due to an unfavourable inclination or anteversion of the cup (acetabulum), which can be compensated by the asymmetric inlay
- Instabilities in the hip joint endoprosthesis due to unfavourable inclination or anteversion of the cup, which can be compensated by the asymmetric inlay

REVISIO® PE-Inlay Constrained

- Recurrent dislocation of the hip joint endoprosthesis due to soft tissue weakness
- Failed alternative methods of treatment for hip joint instability

- Primary or retreatment of patients with a high risk of hip joint dislocation where an alternative therapy does not appear promising, e.g. due to
 - Abductor insufficiency
 - Inadequate bone and/or soft tissue conditions
 - Cognitive or neuromuscular abnormalities
 - Instabilities with unclear aetiology
 - Multidirectional intraoperative instability
- Older patients with a high risk of hip joint dislocation and reduced expectations regarding motion

The primary use of the asymmetric inlay is recommended in exceptional cases where there is a corresponding risk of dislocation but where repositioning of the cup does not appear feasible for stability reasons.

6. Contraindications

6.1 General contraindications

- Acute or chronic, local or systemic, unmanageable infection
- Deficient bone quality or bone conditions that endanger a stable anchorage of the prosthesis
- Vascular, muscle and nerve diseases that endanger the affected area
- Any accompanying disease that can impair bone regeneration or the implant's function (e.g. severe metabolic disorders)
- High physical strain (e.g. competitive sports, heavy physical labour)
- Incomplete skeletal growth
- Allergies to materials of implant components
- Extreme obesity or severe weight gain
- Pregnancy and breastfeeding
- Lack of compliance
- Taking high doses of cortisone or cytostatic drugs
- Epilepsy or other circumstances which could lead to accidents with an increased risk of fractures

Contraindications also include those that generally exclude surgery, such as coagulation disorders (also due to anticoagulants that are indispensable), untreatable high blood pressure, recent heart attack, cardiovascular diseases, severe cardiac insufficiency, intolerance to anaesthesia or

poor general condition of the patient. To date, no reliable scientific findings are available for use during pregnancy and lactation. For this reason, the implant should not be used in pregnant and breastfeeding women. Ongoing chemotherapy or the aggressive treatment of rheumatic diseases require special clarification.

6.2 Specific contraindications

REVISIO® S Acetabular Cup

- Insufficient peripheral anchorage options
- Missing acetabular cup rim (e.g. due to major acetabular cup defects, so that there are no sufficient anchorage options)

REVISIO® M Acetabular Cup Standard

- Insufficient peripheral anchorage options at the Os ilium and Os ischii
- Defects more pronounced than type IIb according to Paprosky

REVISIO® M Acetabular Cup with flange

- Insufficient peripheral anchorage options at the Os ilium and Os ischii
- Defects more extensive than type IIIb according to Paprosky

REVISIO® PE-Inlays Standard

- Increased risk of postsurgical dislocation

REVISIO® PE-Inlays Asymmetric 10°/20° and Constrained

- Primary use in patients without higher risk of dislocation

Contraindications for cementless fixation

- Osteoporotic bone
- Osteodystrophy
- Osteopenia
- Metabolic bone disorders

7. Potential adverse effects

The expected life expectancy of an implant can be impaired in young, very tall and obese patients and patients who are physically very active. Wrong selection or maladjustment of the

prosthesis, improper fixation, positioning or use despite of contraindication, can lead to early failure due to loosening, fracture, or wear.

Typical adverse effects of artificial joints may include:

- Aseptic implant loosening
- Periprosthetic infection
- Chronic pain
- Luxation
- Limited implant survival due to obesity
- Changes of the bone tissue (atrophy, necrosis)
- Hypersensitivity to alloy components
- Limited joint function and mobility

8. Patient information

Serial or batch number as well as item number of the implanted prosthesis should be documented in the patient's record and the implant pass the patient receives after surgery. Labels are included in the sterile packaged implants.

The patient must be informed about potential benefits and risks of the procedure and the expected effects resulting from individual physiological conditions, e.g. bone quality.

The patient must be informed about the proper handling of the implant and avoidance of harmful activities that might lead to implant loosening or fatigue as well as bone fractures. The patient should also be informed that the degree of motion prescribed during the informed consent discussion must be observed. Moreover, the patient should be informed that sports activities are possible to a limited extend and after consultation and individual risk assessment by the attending physician.

Patient information should also include the possible need for partial or total revision surgery, e.g. due to implant loosening as a consequence of mechanical wear or infection.

The patient must be informed that the body weight plays a decisive role in joint replacement and that heavy loads should be avoided.

However, the patient must be informed that regular physiotherapy should be carried out in accordance with the individual instructions of the attending physician.

The patient must be informed that diagnostic or follow-up magnetic resonance imaging (MRI) can cause interactions with the implant.

The patient should be informed that the use of cups or asymmetrical cup inlays can result in a limited range of motion (ROM).

Finally, the patient should be advised in which circumstances the attending physician should be contacted.

9. Implant materials

The following implant materials are used:

- Cobalt alloy (CoCrMo) manufactured according to ISO 5832-4 or ISO 5832-12, respectively (REVISIO® S and M Acetabular Cups)
- UHMWPE – Ultra-high molecular weight polyethylene according to ISO 5834-1, -2 (Inlays)
- Titanium-Niobium [TiNb] coating of CoCrMo components for patients suffering from allergic reactions to nickel, cobalt, chromium, molybdenum
- Titanium alloy (TiAl6V4), produced according to ISO 5832-3 (screws, iliac peg)

Note:

Please follow the instructions regarding the mechanical compatibility of the implant materials and consider potential allergic reactions to one or multiple components of these implant materials (see Contraindications).

10. Clinical follow-up

Patients with an artificial joint require continuous follow-up by the surgeon or another competent physician. In the event of an unforeseen early implant loosening, material weakening or similar signs of potential complications, appropriate measures can be taken in time. If a revision procedure should become necessary, the chances for success are significantly higher in case of early identification, especially if the patient's bone quality at the implant-bone interface has not yet deteriorated significantly.

The patient should be advised to inform the surgeon or attending physician about any changes relating to the artificial joint.

Regular clinical follow-up examinations are recommended. If this is not possible, the surgeon or attending physician should be provided with a control X-ray for assessment at least once a year.

AQ Solutions GmbH provides appropriate instruments for revision procedures. Please refer to the specific instructions in our user information and surgical techniques. Upon request, our product specialists will provide support during revision surgery at your convenience.

Caution:

In case of any queries relating revision procedures, please contact the implant manufacturer.

Caution:

Exposure to strong magnetic or electromagnetic fields commonly used in magnetic resonance imaging (MRI) can have adverse effects on the implant behavior.

11. Compatibility

Caution:

We guarantee the compatibility of our products only in connection with our own products and the products from other manufacturers we have approved for the combination. The combination of AQ Solutions implants with not approved products from other manufacturers is excluded for reasons of product safety and product liability. Please refer to the note on the compatibility of products from AQ Solutions GmbH in our OP-Manual.

For the REVISIO® S and M Acetabular Cups, only locking, reducing and cover screws supplied by AQ Solutions GmbH must be used.

For the REVISIO® M Acetabular Cup, only the iliac peg including locking screw supplied by AQ Solutions GmbH must be used.

12. Instrument care – special precautions

Caution:

Proper handling, decontamination, including pre-rinsing, washing, rinsing and sterilisation, storage, and use of instruments are decisive for a long and adequate durability of all surgical instruments. However, even with correct use, maintenance and care, life expectancy of the instruments is limited. This applies, in particular, to cutting instruments and force-controlled instruments. These instruments and their components are often exposed to high loads and/or impaction forces. Defects can occur under such conditions, especially if the instrument is corroded, damaged, notched, or scratched. Please replace any instruments that have become unusable. Examine all instruments carefully before using them. Never use an instrument that is damaged or worn, or has an illegible laser engraving or a cutting instrument with blunt cutting edges.

Note:

Please follow the general instructions on the maintenance and care of surgical instruments.

13. Symbols

Appendix 2 contains the symbols used on the package labels.

Bezeichnung • Designation	Artikelnummer • Item number		
REVISIO® S PFANNE REVISIO® S ACETABULAR CUP	22010344 22025350 22040356 22060364	22015346 22030352 22045358 22070368	22020348 22035354 22050360
REVISIO® M PFANNE STANDARD REVISIO® M ACETABULAR CUP STANDARD	20007044 20007056 20007068	20007048 20007060 20007072	20007052 20007064 20007076
REVISIO® M PFANNE MIT LASCHE REVISIO® M ACETABULAR CUP WITH FLANGE	20009052 20009064 20009076 20009160 20009172	20009056 20009068 20009152 20009164 20009176	20009060 20009072 20009156 20009168
REVISIO® S/M VERSCHLUSSSCHRAUBE REVISIO® S/M COVER SCREW	20000108	20000408	
REVISIO® VERSCHLUSSSCHRAUBE PFANNENPOL REVISIO® POLE COVER SCREW	20000110		
REVISIO® VERSCHLUSSSCHRAUBE REVISIO® COVER SCREW	20011014		
REVISIO® REDUZIERSCHRAUBE REVISIO® REDUCING SCREW	20012014		
REVISIO® SICHERUNGSSCHRAUBE DARMBEINZAPFEN REVISIO® LOCKING SCREW	20013014		
REVISIO® DARMBEINZAPFEN REVISIO® FIXATION PEG	20010030	20010050	20010070

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