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AQ SOLUTIONS GMBH GUIDELINES FOR TREATING REUSABLE SURGICAL INSTRUMENTS OF AQ SOLUTIONS GMBH

according to DIN EN ISO 17664

This guide must be read carefully before using instruments made by AQ Solutions GmbH!

1. GENERAL PRINCIPLES

These guidelines apply for all AQ Solutions GmbH instruments. Furthermore, we recommend classifying these instruments into Critical A and / or Critical B for the respective instruments according to KRINKO.

All instruments must be cleaned, disinfected, and sterilised prior to each application; this is also required prior to the first use after delivery of the non-sterile instruments (cleaning and disinfection after removal of the protective packaging, sterilisation after packaging). Effective cleaning and disinfection are mandatory for an effective sterilisation of the instruments.

You are responsible for the sterility of the instruments. Therefore, please ensure

- that only procedures validated for the devices and products are used for cleaning, disinfection, and sterilisation
- that all used devices (WD, steriliser) are maintained and reviewed regularly
- the validated parameters comply with for every cycle

Make sure to avoid heavy contamination of the instrument tray during application; otherwise it is necessary to clean and disinfect the instruments and the instrument tray.

Moreover, please pay attention to the national legal provisions as well as to the hygienic instructions of the respective practice or hospital. This applies particularly to the different guidelines regarding the inactivation of prions.

Instrument care – special precautions

Proper handling, decontamination, including pre-rinsing, washing, rinsing and sterilisation, storage, and use of instruments are decisive for a long and adequate durability of all surgical instruments. However, even with correct use, maintenance and care, life expectancy of the instruments is limited. This applies, in particular, to cutting instruments and force-controlled instruments. These instruments and their components are often exposed to high loads and/or impaction forces. Defects can occur under such conditions, especially if the instrument is corroded, damaged, notched, or scratched. Please replace any instruments that have become unusable. Examine all instruments carefully before using them. Never use an instrument that is damaged or worn, or has an illegible laser engraving or a cutting instrument with blunt cutting edges.

Please note: For any additional or deviating provisions please refer to Appendix 2 'Specific instructions'.

2. CLEANING AND DISINFECTION

Basic principles

Where appropriate, automated procedures (WD, washer-disinfector) should be used for cleaning and disinfection of the instruments. Due to the significantly lower efficiency and reproducibility of manual procedures – even in combination with an ultrasonic bath – they should only be used if an automated procedure is not available.

Pre-treatment

Please remove heavy dirt directly after application (within a maximum of 2 hours).

- Procedure:
1. Disassemble instruments as far as possible.
 2. Rinse instruments at least 1 min under running water (temperature < 35 °C/95 °F)
Rinse all lumens at least three times with a syringe (minimum volume of 50ml).
Move movable parts back and forth at least three to five times during pre-rinsing.
 3. Immerse disassembled instrument parts in the pre-cleaning solution (ultrasonic bath with deactivated ultrasonics) using a small components basket for the specified soaking time ensuring that the instruments are sufficiently covered (see Appendix 2 'Specific instructions'). Pay attention that the components do not come into contact with each other. Support the cleaning process by thoroughly brushing of all outer and inner surfaces with a soft brush (at beginning of soaking, please open the fine mesh wire basket for this and close afterwards). Move movable parts back and forth at least three to five times during pre-cleaning). Rinse all lumens of the instruments at least three times at the beginning and at the end of the soaking time with a syringe (minimum volume of 50ml).
 4. Activate ultrasonics for an additional soaking time (at least 15 min).
 5. Remove instruments and rinse them at least three times thoroughly (at least 1 min) with water.
Rinse all lumens of the instruments at least three times at the beginning and at the end of the soaking time with a syringe (minimum volume of 50ml).
Move movable parts back and forth at least three to five times during post-rinsing.

 **Caution: Some instruments require additional or deviating provisions; please refer to Appendix 2 'Specific instructions'.**

For the selection of the detergent, the following aspects should be considered¹:

- Suitability for the cleaning of instruments made of metal or plastic
- Suitability for ultrasonic cleaning (no foam development)
- Compatibility with the instruments (see chapter 'Material durability')

The manufacturer's instructions regarding concentration, temperature, soaking time, and post-rinsing must be observed. Please use only freshly prepared solutions and sterile or low-germ water (max. 10 germs/ml), as well as low-endotoxin water (max. 0.25 endotoxin units/ml), e.g. purified/highly purified water, a soft and clean lint-free cloth and/or filtered air for drying, respectively.

3. AUTOMATED CLEANING/DISINFECTION (WD, WASHER-DISINFECTOR)

For the selection of the WD, the following aspects should be considered:

- proven efficacy (e.g. CE marking according to EN ISO 15883 or DGHM approval)
- possibility to use an approved program for thermal disinfection (A_0 value $\geq 3,000$ or – in case of older devices – at least 5 min at 90 °C/194 °F; in case of chemical disinfection, there is a risk of disinfectant residues remaining on the instruments)
- suitability of the program for instruments and sufficient rinsing steps in the program
- post-rinsing only with sterile or low-germ water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), e.g. purified/highly purified water
- air used for drying is filtered (oil-free, low contamination with microorganisms and particles)
- regular maintenance and review/calibration of the WD

For the selection of the detergent system, the following aspects should be considered:

- suitability for the cleaning of instruments made of metal or plastic
- if thermal disinfection is not used, a suitable disinfectant with proven efficacy (e.g. VAH/DGHM approval or CE marking) and compatible with the respective detergent is used as an alternative
- compatibility of the used detergents with the instruments (see chapter Material durability¹)

The manufacturer's instructions regarding concentration, temperature and soaking time, as well as post-rinsing must be observed.

Procedure:

1. Disassemble instruments as far as possible.
2. Place disassembled parts in the WD using a small components basket, where applicable (see Appendix 2 'Specific instructions'). Pay attention that the components do not come into contact with each other. Connect all lumens to the rinsing port of the WD.
3. Start the program.
4. Disconnect and remove instruments from the WD after the end of the program.
5. After subsequent drying in a clean area (where applicable) review and package the instruments immediately after removal (see chapters 'Review', 'Maintenance', and 'Packaging').

Suitability of the instruments for effective automated cleaning and disinfection was proven by an independent, officially accredited and recognised (§ 15 (5) MPG) test laboratory using WD G 7836 CD, Miele & Cie. GmbH & Co., Gütersloh, (thermal disinfection) and Neodisher® medizym (Dr Weigert GmbH & Co. KG, Hamburg) using the above mentioned procedure.

4. MANUAL CLEANING AND DISINFECTION

For the selection of the detergent and disinfectant, the following aspects should be considered:

- suitability for the cleaning and disinfection of instruments made of metal or plastic
- when using an ultrasonic bath: suitability of the detergent for ultrasonic cleaning (no foam development)
- use of a disinfectant with proven efficacy (e.g. VAH/DGHM approval or CE marking) compatible with the used detergent
- compatibility of the chemicals with the instruments (see chapter 'Material durability'¹)

¹ If a cleaning detergent and disinfectant is used (e.g. for occupational health reasons) please consider, that it should be free of aldehydes (to avoid adhesion of blood particles) with proven efficiency (e.g. VAH/DGHM approval or CE marking) suitable for the disinfection of instruments made of metal or plastic, and compatible with the instruments (see chapter 'Material durability'). Please note that a disinfectant is only used during pre-treatment for the employee's personal safety, but cannot replace the disinfection step after cleaning.

The manufacturer's instructions regarding concentration, temperature and soaking time for detergent and disinfectant as well as post-rinsing must be observed. Please use only freshly prepared solutions and sterile or low-germ water (max. 10 germs/ml) as well as low-endotoxin water (max. 0.25 endotoxin units/ml), e.g. purified/highly purified water, a soft and clean lint-free cloth and/or filtered air (oil-free, low contamination regarding germs and particles) for drying, respectively.

 **Notice: In general, all products in Appendix 1 should not be cleaned manually due to the lower log reduction.**

Procedure: Cleaning

1. Disassemble instruments as far as possible.
2. Immerse disassembled instrument parts in the cleaning solution for the specified soaking time using a small components basket, where applicable (see Appendix 2 'Specific instructions'). Pay attention that the components do not come in contact with each other. Support the cleaning process by thoroughly brushing of all outer and inner surfaces with a soft brush (at beginning of soaking, please open the fine mesh wire basket for this and close afterwards). Move movable parts back and forth at least three to five times during cleaning. Rinse all lumens of the instruments at least three times at the beginning and at the end of the soaking time with a syringe (minimum volume of 50ml).
3. Activate ultrasonics for an additional soaking time (at least 15 min).
4. Remove instruments and rinse them at least three times thoroughly (at least 1 min) with water. Move movable parts back and forth at least three to five times during post-rinsing. Rinse all lumens at least five times at the beginning and at the end of the soaking time with a single-use syringe (minimum volume of 50ml).
5. Review instruments (see chapters 'Review' and 'Maintenance').

Disinfection

6. Immerse the disassembled instrument parts in the disinfectant solution for the specified soaking time so that the instruments are sufficiently covered. Pay attention that the components do not come in contact with each other. Move movable parts back and forth at least three to five times during disinfection. Rinse all lumens at least five times at the beginning and at the end of the soaking time with a single-use syringe (minimum volume of 50ml)
7. Remove instruments and rinse them at least five times thoroughly (at least 1 min) with water. Rinse all lumens at least five times at the beginning and at the end of the soaking time with a single-use syringe (minimum volume of 50ml). Move movable parts back and forth at least three to five times during post-rinsing.
8. Dry (with filtered compressed air) and package instruments immediately after removal (see chapter 'Packaging', if necessary, after additional subsequent drying at a clean place).

Suitability of the instruments for an effective cleaning and disinfection was proven by an independent, officially accredited and recognised (§ 15 (5) MPG) test laboratory using the pre-cleaning and detergent Cidezyme/Enzol and disinfectant Cidex® OPA (Johnson & Johnson GmbH, Norderstedt) and the above mentioned procedure.

5. REVIEW

After cleaning or cleaning/disinfection, review all instruments for potential corrosion, damaged surfaces, chipping, contamination, and discoloration. Do not use damaged instruments (for the maximum of reuse cycles refer to chapter 'Reusability'). Any instruments that still show signs of contamination must be cleaned and disinfected again.

6. MAINTENANCE/ASSEMBLY

Assembly is only permitted after cleaning/disinfection.
Modular instrument parts should be assembled immediately prior to use.
Instrument oil or grease shall, in principle, not be used.

7. PACKAGING

Place the cleaned and disinfected instruments into the corresponding sterilisation trays and sterilisation containers complying with the following requirements (material/process):

- EN ISO 11607 (in accordance with EN 868-8)
- suitable for steam sterilisation (temperature resistant up to at least 138 °C (280 °F), sufficient steam permeability)
- adequate protection of instruments and sterilisation packaging against mechanical damage
- regular maintenance according to the instructions of the manufacturer (sterilisation container)

8. STERILISATION

Only listed sterilisation procedures should be used for sterilisation; other sterilisation procedures are not permitted.

Steam sterilisation

- fractionated vacuum/dynamic air removal procedure^{2,3} (with adequate product drying⁴)
- steam steriliser according to EN 285/EN 13060
- validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance qualification (PQ))
- maximum sterilisation temperature of 134 °C (273 °F; plus tolerance according to EN ISO 17665)
- sterilisation time (exposure time at the sterilisation temperature):

Area	Fractionated vacuum/dynamic air removal	Gravity procedure
Germany	at least 5 min at 134 °C ⁴	not admitted
other countries	at least 5 min ⁵ at 132 °C (270 °F) / 134 °C (273 °F), minimum drying time of 20 min ⁴	not admitted

Suitability of the instruments for an effective steam sterilisation was proven by an independent, officially accredited and recognised (§ 15 (5) MPG) test laboratory using the steam steriliser HST 6x6x6 (Zirbus technology GmbH, Bad Grund) and fractionated vacuum/dynamic air removal procedure. Therefore, normal conditions at the clinic and the specified procedure were taken into consideration.

Flash sterilisation is prohibited.

Do not use dry heat sterilisation, radiation sterilisation, formaldehyde and ethylene oxide sterilisation, as well as plasma sterilisation.

² at least three vacuum steps

³ The less effective gravity procedure shall only be used if the fractionated vacuum procedure is not possible; it requires significantly longer sterilisation times as well as device-, procedure-, parameter-, and product-specific validation that are the sole responsibility of the user.

⁴ The required drying time depends on parameters that are the sole responsibility of the user (load configuration and density, steriliser conditions, ...) and are thus to be determined by the user. Nevertheless, drying time should not be below 20 min.

⁵ respectively 18 min (inactivation of prions)

9. STORAGE

After sterilisation, instruments are stored at a dry and dust-free place in the sterilisation packaging. The provisions of DIN EN 11607 shall apply.

10. MATERIAL DURABILITY

Make sure that the cleaning or disinfection detergents do not contain the following ingredients:

- organic, mineral, and oxidising acids (minimum permitted pH-value of 5.5)
- organic solvents (e.g. alcohol, ether, ketone, petroleum)
- oxidising agents (e.g. peroxide)
- Halogens (e.g. chlorine, iodine, bromine)
- aromatic, halogenated hydrocarbons
- Ammonia

Please do not clean instruments and sterilisation trays with metal brushes or steel wool.

Because of the chemical resistance of the material Ferrozell, using an alkaline cleaning and disinfection agent on it may cause undesirable changes to the material.

It is recommended to use a neutral, enzymatic cleaning agent in order to retain the material properties of components made of Ferrozell. If alkaline cleaning agents are used, however, and this causes changes to the material, the instruments in question have to be discarded (see point 5 and 11).

Please do not expose instruments and sterilisation trays to temperatures exceeding 138 °C (280 °F).

11. REUSABILITY

All instruments can be reused, if they are handled carefully and are movable, undamaged, and clean. Any further use or the use of damaged and dirty instruments shall be the sole responsibility of the user.

For instruments containing ferrozell (wood-like structure), care must be taken to ensure that the surface is not damaged (cracks, gaps between ferrozell and metal). In case of visible damage of the surface, the instruments must be replaced. No liability in case of disregard.

TRADEMARKS AND DISCLAIMERS

This document is intended for medical professionals only. The surgeon must always use his/her own professional clinical judgment in deciding whether or not to use a particular product in the treatment of a patient. AQ Solutions does not provide medical advice and recommends that surgeons be trained in the use of a particular product before using that product in a procedure. Surgeons must always refer to the package insert, product label and/or instructions for use, including cleaning and sterilization instructions, if applicable, before using an AQ Solutions product. Please contact your AQ Solutions representative if you have any questions regarding the use and availability of any of our products. AQ Solutions is the owner of the following trademarks: REVISIO®.

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Artikel-Nr.	Bezeichnung	Designation
02000025	Hüftpfannen-Setzinstrument, gebogen	Cup impactor, curved
02004903	Tiefenmesslehre gebogen	Depth gauge curved
02005025	Kardanschraubendreher SW3,5 - Lang, Dreiecksgriff	Flexible screwdriver SW3,5 - Long, triangular handle
02003112	Schraubendreher - Reduzierschraube	Screwdriver - reducing screw
02005100	Hüftkopf Eindrückinstrument PE-Inlay constrained	Cup impressor PE-Inlay constrained
02000100	T-Handgriff AO	T-handle AO
02004052 - 02004176	Probepfanne REVISIO® M mit Lasche Links/Rechts 52 mm – 76 mm	Trial cup REVISIO® M with flange left/right 52 mm – 76 mm
02904544	Tibia-Höhenfühler medial/lateral	Gauge tibial resection guide medial/lateral
02900161	Nagelrevisionszange	Pin puller
02904653	Halter für Tibia-Größenlehre	Handle for tibial sizer
02904594	Femur-Messlehre	Femoral gauge
03000098	Femurfühler anterior S	Femoral probe anterior S
03000099	Femurfühler anterior M	Femoral probe anterior M
03000100	Femurfühler anterior L	Femoral probe anterior L
02905466	Femur-Sägelehre 4 in 1 S	Femoral cutting guide 4 in 1 S
02905467	Femur-Sägelehre 4 in 1 M	Femoral cutting guide 4 in 1 M
02905468	Femur-Sägelehre 4 in 1 L	Femoral cutting guide 4 in 1 L
04000201	Femur-Sägelehre Adapter 4 in 1 S	Femoral cutting guide adapter 4 in 1 S
04000202	Femur-Sägelehre Adapter 4 in 1 M	Femoral cutting guide adapter 4 in 1 M
04000203	Femur-Sägelehre Adapter 4 in 1 L	Femoral cutting guide adapter 4 in 1 L
02904519	Femurkasten-Sägelehre S	Femoral box cutting guide S
02904520	Femurkasten-Sägelehre M	Femoral box cutting guide M
02904521	Femurkasten-Sägelehre L	Femoral box cutting guide L
02905816	Ein-/Ausschläger Femurkasten-Sägelehre	Impactor/extractor for femoral box cutting guide
02905765	T-Handgriff Formraspel	T-handle for rasp
00200734	Raspelhandgriff Endofemoral	Rasp handle endofemoral
03000300	Raspelhandgriff Endofemoral	Rasp handle endofemoral
08000101	Gleithammer T-Griff	Slide hammer T-Handle

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Geometry aspect	Special/Additional procedure in case of			
	Pre-treatment	Manual cleaning/ disinfection	Automated cleaning/ disinfection	Sterilisation
Small parts	Use of small components-basket with lid	Use of small components basket with lid	Use of small components basket with lid	—
Long cannulation / lumen	Always rinse with single-use syringe	Always rinsing with single-use syringe	<ul style="list-style-type: none"> • mounted to a rinse lance • Connection to a rinse tube 	—
Rasp	Mandatory visual control at the end of pre-cleaning, in case of remnants repetition of pre-cleaning and discard instrument	Mandatory visual control at the end of cleaning, in case of remnants repetition of cleaning and discard instrument	Mandatory visual control at the end of cleaning/disinfection, in case of remnants repetition of cleaning/disinfection and discard instrument	—
Blind hole	Always brush out and back-flush with single-use syringe and cannula	Always brush out and back-flush with single-use syringe and cannula	Always with opening downwards or sideways	—
Screw sleeve	Always with unfastened screw sleeve	Always with unfastened screw sleeve	Always with unfastened screw sleeve	—
O-ring	Disassemble	Disassemble	Disassemble	Assemble
Joint	Always with open joint	Always with open joint	Always with open joint	—
Multi-part instruments	Disassemble so that all surfaces are accessible for cleaning and disinfection; in case of complex instruments observe the specific assembly instructions	Disassemble so that all surfaces are accessible for cleaning and disinfection; in case of complex instruments observe the specific assembly instructions	Disassemble so that all surfaces are accessible for cleaning and disinfection; in case of complex instruments observe the specific assembly instructions	Disassemble
Slot	Brush with a small pipe brush	Brush with a small pipe brush	—	—
Silicon strap	Disassemble	Disassemble	Disassemble	Disassemble
Locking screw	Disassemble	Disassemble	Disassemble	Disassemble
Cardan screwdriver	Cardan screwdrivers must not be cleaned in an ultrasonic bath.			—
T-Handle AO	—			Apply water-soluble instrument lubricant to moving parts

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